

For Research Use Only

Mouse Il13ra1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87816-1-PBS



Basic Information

Catalog Number: 87816-1-PBS	GenBank Accession Number: NM_133990.5	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 16164	CloneNo.: 253331D8
Isotype: IgG	UNIPROT ID: O09030	
Immunogen Catalog Number: EG6217	Full Name: interleukin 13 receptor, alpha 1	
	Calculated MW: 48 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IL-13Ra1/CD213a1 concentration was determined to be 21.9 ng/mL in mouse liver tissue extract based on a 8.4 mg/mL extract load, mouse skeletal muscle tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03221-1, MOUSE IL-13Ra1/CD213a1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87816-2-PBS. Detection antibody: 87816-1-PBS. Standard: Eg6217. Range: 0.312-40 ng/mL.