

For Research Use Only

MU154079-4 (ATXN3) Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87817-2-PBS



Basic Information

Catalog Number:

87817-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG41832

GenBank Accession Number:

GeneID (NCBI):

4287

Full Name:

ataxin 3

Purification Method:

Protein A purification

CloneNo.:

253467E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Ataxin-3/ATXN3 concentration was determined to be 567.7 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load, 353.4 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load, 132.2 ng/mL in U-251 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03223-1, Ataxin-3/ATXN3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87817-3-PBS. Detection antibody: 87817-2-PBS. Standard: Ag41832. Range: 1.562-100 ng/mL.