

For Research Use Only

FCRL2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87834-2-PBS



Basic Information

Catalog Number: 87834-2-PBS	GenBank Accession Number: NM_030764.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 79368	CloneNo.: 253006C3
Isotype: IgG	UNIPROT ID: Q96LA5-1	
Immunogen Catalog Number: EG4773	Full Name: Fc receptor-like 2	
	Calculated MW: 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FCRL2 concentration was determined to be 14.6 ng/mL in Daudi cell extract based on a 1.2 mg/mL extract load, 9.5 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03222-1, FCRL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87834-2-PBS. Detection antibody: 87834-1-PBS. Standard: Eg4773. Range: 0.625-80 ng/mL.