

For Research Use Only

FCRL2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87834-3-PBS



Basic Information

Catalog Number:

87834-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4773

GenBank Accession Number:

NM_030764.3

GeneID (NCBI):

79368

UNIPROT ID:

Q96LA5-1

Full Name:

Fc receptor-like 2

Calculated MW:

56 kDa

Purification Method:

Protein A purification

CloneNo.:

253006C2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



The mean FCRL2 concentration was determined to be 6.75 ng/mL in Daudi cell extract based on a 1.2 mg/mL extract load and 5.41 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load.



Plasma of six individual healthy human donors was measured. The human FCRL2 concentration of detected samples was determined to be 1,294.6 pg/mL with a range of 980.1-1,641.9 pg/mL.



Sandwich ELISA standard curve of MP03222-2, Human FCRL2 Recombinant Matched Antibody Pair - PBS only. 87834-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4773. 87834-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL