

For Research Use Only

BMI1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87836-3-PBS



Basic Information

Catalog Number:

87836-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1286

GenBank Accession Number:

BC011652

GeneID (NCBI):

648

UNIPROT ID:

P35226

Full Name:

BMI1 polycomb ring finger oncogene

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

253339D8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BMI 1 concentration was determined to be 7.4 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load, 4.4 ng/mL in NIH/3T3 cell extract based on a 1.5 mg/mL extract load, 5.1 ng/mL in U-251 cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP03232-1, BMI1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87836-3-PBS. Detection antibody: 87836-2-PBS. Standard: Ag1286. Range: 1.25-80 ng/mL.