

For Research Use Only

TNFRSF6B Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87840-2-PBS



Basic Information

Catalog Number:

87840-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3110

GenBank Accession Number:

NM_003823.3

GeneID (NCBI):

8771

UNIPROT ID:

O95407

Full Name:

tumor necrosis factor receptor
superfamily, member 6b, decoy

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

253457D10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



SW 480 cells were cultured in DMEM supplemented with 10% fetal bovine serum, 2.5 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. An aliquot of the cell culture supernate was removed, assayed for human DcR3/TNFRSF6B, and measured 29.05 ng/mL.



Plasma of five individual healthy human donors was measured. The human DcR3/TNFRSF6B concentration of detected samples was determined to be 2.47 ng/mL with a range of 0.38 - 9.39 ng/mL.



Sandwich ELISA standard curve of MP03265-1, Human DcR3/TNFRSF6B Recombinant Matched Antibody Pair - PBS only. 87840-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3110. 87840-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL