

For Research Use Only

CWC15 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87890-2-PBS



Basic Information

Catalog Number: 87890-2-PBS	GenBank Accession Number: BC040946	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 51503	CloneNo.: 260028G9
Isotype: IgG	UNIPROT ID: Q9P013	
Immunogen Catalog Number: AG18116	Full Name: CWC15 spliceosome-associated protein homolog (S. cerevisiae)	
	Calculated MW: 229 aa, 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CWC15 concentration was determined to be 7.5 ng/mL in Ramos cell extract based on a 1.2 mg/mL extract load, 5.7 ng/mL in PC-3 cell extract based on a 1.2 mg/mL extract load, 7.9 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03260-1, CWC15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87890-1-PBS. Detection antibody: 87890-2-PBS. Standard: Ag18116. Range: 0.156-20 ng/mL.