

For Research Use Only

CCT3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87920-1-PBS



Basic Information

Catalog Number:

87920-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0843

GenBank Accession Number:

BC006501

GeneID (NCBI):

7203

UNIPROT ID:

P49368

Full Name:

chaperonin containing TCP1, subunit 3 (gamma)

Calculated MW:

60 kDa

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

260042D9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

CCT3 protein (full name: chaperonin containing TCP1 subunit 3, also known as the gamma subunit of TCP-1 protein) is an important member of the TCP-1 chaperonin family. As one of the eight essential subunits of the CCT/TRiC complex, its core function is to act as a molecular chaperone that assists in the correct folding of other cellular proteins under ATP hydrolysis. It is particularly involved in the folding process of key cytoskeletal components such as actin and tubulin, while also participating in important biological activities including ciliogenesis and telomere maintenance.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87920-1-RR (CCT3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87920-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87920-1-RR against Human CCT3 were performed. The affinity constant is below 1 pM.