

For Research Use Only

STAC3 Recombinant monoclonal antibody

Catalog Number: 87924-1-RR



Basic Information

Catalog Number:

87924-1-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5527

GenBank Accession Number:

NM_145064.3

GeneID (NCBI):

246329

UNIPROT ID:

Q96MF2-1

Full Name:

SH3 and cysteine rich domain 3

Calculated MW:

42 kDa

Observed MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

260138A8

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, pig

Positive Controls:

WB : human skeletal muscle tissue, mouse skeletal muscle tissue, pig skeletal muscle tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87924-1-RR (STAC3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.