

For Research Use Only

FAM171A1 Recombinant monoclonal antibody

Catalog Number:87948-1-RR



Basic Information

Catalog Number: 87948-1-RR	GenBank Accession Number: NM_001081161.1	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 269233	CloneNo.: 260116C2
Isotype: IgG	UNIPROT ID: A2ATK9	Recommended Dilutions: WB: 1:5000-1:50000
Immunogen Catalog Number: EG8071	Full Name: family with sequence similarity 171, member A1	
	Calculated MW: 98 kDa	
	Observed MW: 100-120 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse brain tissue, rat brain tissue
Species Specificity: mouse, rat	

Background Information

FAM171A1 is a transmembrane glycoprotein that regulates actin dynamics and promotes tumor invasiveness in several cancers while also showing genetic associations with neurodegenerative diseases.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87948-1-RR (FAM171A1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.



mouse brain tissue were subjected to SDS PAGE followed by western blot with 87948-1-RR (FAM171A1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.