

For Research Use Only

PIK3C3 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87952-2-PBS



Basic Information

Catalog Number:

87952-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3110

GenBank Accession Number:

BC033004

GeneID (NCBI):

5289

UNIPROT ID:

Q8NEB9

Full Name:

phosphoinositide-3-kinase, class 3

Calculated MW:

887 aa, 100 kDa

Purification Method:

Protein A purification

CloneNo.:

260118A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean VPS34 concentration was determined to be 54.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 15.6 ng/mL in PC-3 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03292-1, VPS34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87952-1-PBS. Detection antibody: 87952-2-PBS. Standard: Ag3110. Range: 0.312-40 ng/mL.