

For Research Use Only

PSG7 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87954-1-PBS



Basic Information

Catalog Number:

87954-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG8784

GenBank Accession Number:

NM_002783.3

GeneID (NCBI):

5676

UNIPROT ID:

Q13046

Full Name:

pregnancy specific beta-1-glycoprotein 7

Calculated MW:

47 kDa

Purification Method:

Protein A purification

CloneNo.:

260017B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PSG7 concentration was determined to be 3262.7 ng/mL in maternal plasma.



Cytometric bead array standard curve of MP03280-1, PSG7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87954-2-PBS. Detection antibody: 87954-1-PBS. Standard: Eg8784. Range: 0.156-20 ng/mL.