

For Research Use Only

PSG5 Recombinant monoclonal antibody

Catalog Number: 87956-1-RR



Basic Information

Catalog Number: 87956-1-RR	GenBank Accession Number: NM_002781.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5673	CloneNo.: 260012B6
Isotype: IgG	UNIPROT ID: Q15238	Recommended Dilutions: WB: 1:2000-1:10000
Immunogen Catalog Number: EG8786	Full Name: pregnancy specific beta-1-glycoprotein 5	
	Calculated MW: 38 kDa	
	Observed MW: 35-45 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human placenta tissue,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87956-1-RR (PSG5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.