

For Research Use Only

MFF Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87972-2-PBS



Basic Information

Catalog Number:

87972-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9990

GenBank Accession Number:

BC000797

GeneID (NCBI):

56947

UNIPROT ID:

Q9GZY8

Full Name:

mitochondrial fission factor

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

260043B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MFF concentration was determined to be 75.4 ng/mL in Hela cell extract based on a 1.2 mg/mL extract load, 0.8 ng/mL in BxPC-3 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03305-1, MFF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87972-3-PBS. Detection antibody: 87972-2-PBS. Standard: Ag9990. Range: 0.312-40 ng/mL.