

For Research Use Only

# ARHGAP28 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87974-2-PBS



## Basic Information

Catalog Number: 87974-2-PBS	GenBank Accession Number: BC065274	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 79822	CloneNo.: 253589A3
Isotype: IgG	UNIPROT ID: B4DXL2	
Immunogen Catalog Number: AG25293	Full Name: Rho GTPase activating protein 28	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean ARHGAP28 concentration was determined to be 2.0 ng/mL in U2OS cell extract based on a 1.5 mg/mL extract load, 0.3 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03311-1, ARHGAP28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87974-3-PBS. Detection antibody: 87974-2-PBS. Standard: Ag25293. Range: 0.078-20 ng/mL.