

For Research Use Only

ARHGAP28 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87974-3-PBS



Basic Information

Catalog Number: 87974-3-PBS	GenBank Accession Number: BC065274	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 79822	CloneNo.: 253589F5
Isotype: IgG	UNIPROT ID: B4DXL2	
Immunogen Catalog Number: AG25293	Full Name: Rho GTPase activating protein 28	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ARHGAP28 concentration was determined to be 2.0 ng/mL in U2OS cell extract based on a 1.5 mg/mL extract load, 0.3 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03311-1, ARHGAP28 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87974-3-PBS. Detection antibody: 87974-2-PBS. Standard: Ag25293. Range: 0.078-20 ng/mL.