

For Research Use Only

TRPC5 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87978-1-PBS



Basic Information

Catalog Number:

87978-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23271

GenBank Accession Number:

BC137271

GeneID (NCBI):

7224

UNIPROT ID:

Q9UL62

Full Name:

transient receptor potential cation channel, subfamily C, member 5

Calculated MW:

111 kDa

Observed MW:

100-111 kDa

Purification Method:

Protein A purification

CloneNo.:

260014C10

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat, pig

Background Information

Short transient receptor potential channel 5 (TRPC5) is abundantly expressed in hippocampal neurons (PMID: 9553080). TRPC5 channels participate in a variety of physiological and pathophysiological processes, covering vascular smooth muscle activity, endothelial function, adiponectin modulation, oxidative stress, epilepsy, anxiety, pain sensation, as well as cardiac remodeling (PMID: 15538613; 22668831; 16685599; 26069213).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87978-1-RR (TRPC5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87978-1-PBS in a different storage buffer formulation.



IP result of anti-TRPC5 (IP:87978-1-RR, 4ug; Detection:87978-1-RR 1:800) with rat brain tissue lysate 1720 ug. This data was developed using the same antibody clone with 87978-1-PBS in a different storage buffer formulation.