

For Research Use Only

PRPF8 Recombinant monoclonal antibody, PBS Only

Catalog Number:87982-1-PBS



Basic Information

Catalog Number:

87982-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1652

GenBank Accession Number:

BC064370

GeneID (NCBI):

10594

UNIPROT ID:

Q6P2Q9

Full Name:

PRP8 pre-mRNA processing factor 8 homolog (*S. cerevisiae*)

Calculated MW:

274 kDa

Observed MW:

~250 kDa

Purification Method:

Protein A purification

CloneNo.:

260212F12

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87982-1-RR (PRPF8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87982-1-PBS in a different storage buffer formulation.