

For Research Use Only

# PRPF8 Recombinant monoclonal antibody

Catalog Number: 87982-1-RR



## Basic Information

<b>Catalog Number:</b> 87982-1-RR	<b>GenBank Accession Number:</b> BC064370	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 10594	<b>CloneNo.:</b> 260212F12
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> Q6P2Q9	<b>Recommended Dilutions:</b> WB: 1:5000-1:50000
<b>Immunogen Catalog Number:</b> AG1652	<b>Full Name:</b> PRP8 pre-mRNA processing factor 8 homolog ( <i>S. cerevisiae</i> )	
	<b>Calculated MW:</b> 274 kDa	
	<b>Observed MW:</b> ~250 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b>
<b>Species Specificity:</b> human, rat	<b>WB :</b> HeLa cells, HeLa, Jurkat, K-562, A-431, HEK-293T, NIH/3T3, C6 cell

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87982-1-RR (PRPF8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.