

For Research Use Only

Mouse Integrin Alpha 3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88002-1-PBS



Basic Information

Catalog Number:

88002-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4683

GenBank Accession Number:

NM_013565.3

GeneID (NCBI):

16400

UNIPROT ID:

Q62470-1

Full Name:

integrin alpha 3

Calculated MW:

117 kDa

Observed MW:

130-140 kDa

Purification Method:

Protein A purification

CloneNo.:

260079F2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

mouse

Background Information

Integrin Alpha 3 (ITGA3), also named as CD49c, GAPB3, and FRP-2, belongs to the integrin alpha chain family. It is a receptor for fibronectin, laminin, collagen, epiligrin, thrombospondin, and CSPG4. ITGA3 may mediate with LGALS3 the stimulation by CSPG4 of endothelial cell migration. ITGA3 is highly concentrated in epithelial cells, where it strongly adheres to Laminin-5, and Laminin-5-induced rapid adhesion can be blocked by antibodies against the alpha-3 integrin subunit.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 88002-1-RR (Integrin alpha-3/CD49c antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88002-1-PBS in a different storage buffer formulation.



mouse colon tissue were subjected to SDS PAGE followed by western blot with 88002-1-RR (Integrin alpha-3/CD49c antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88002-1-PBS in a different storage buffer formulation.



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 88002-1-RR (Integrin Alpha 3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88002-1-PBS in a different storage buffer formulation.