

For Research Use Only

ERG Recombinant monoclonal antibody, PBS Only

Catalog Number:88009-1-PBS



Basic Information

Catalog Number: 88009-1-PBS	GenBank Accession Number: BC040168	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2078	CloneNo.: 253522H8
Isotype: IgG	UNIPROT ID: P11308	
Immunogen Catalog Number: AG5708	Full Name: v-ets erythroblastosis virus E26 oncogene homolog (avian)	
	Calculated MW: 55 kDa	
	Observed MW: 55 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse

Background Information

ERG (ETS-related gene) is a member of the ETS family of transcription factors, which are downstream effectors of mitogenic signaling transduction pathway. ERG is a key regulator of cell proliferation, differentiation, angiogenesis, inflammation, and apoptosis. ERG is involved in chromosomal translocations, resulting in different fusion gene products, such as TMPSSR2-ERG and NDRG1-ERG in prostate cancer, EWS-ERG in Ewing's sarcoma and FUS-ERG in acute myeloid leukemia.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88009-1-RR (ERG antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88009-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed MOLT-4 cells using ERG antibody (88009-1-RR, Clone: 253522H8) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 88009-1-PBS in a different storage buffer formulation.