

For Research Use Only

# HNRNPH1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88017-2-PBS



## Basic Information

<b>Catalog Number:</b> 88017-2-PBS	<b>GenBank Accession Number:</b> BC001348	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 3187	<b>CloneNo.:</b> 260181A12
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P31943	
<b>Immunogen Catalog Number:</b> AG6506	<b>Full Name:</b> heterogeneous nuclear ribonucleoprotein H1 (H)	
	<b>Calculated MW:</b> 49 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean HNRNPH1 concentration was determined to be 21.9 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 0.3 ng/mL in Y79 cell extract based on a 1.1 mg/mL extract load.



Cytometric bead array standard curve of MP03314-1, HNRNPH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88017-1-PBS. Detection antibody: 88017-2-PBS. Standard: Ag6506. Range: 0.156-20 ng/mL.