

For Research Use Only

Mouse RAMP2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88020-1-PBS



Basic Information

Catalog Number: 88020-1-PBS	GenBank Accession Number: NM_019444.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 54409	CloneNo.: 260084F11
Isotype: IgG	UNIPROT ID: Q9WUP0	
Immunogen Catalog Number: EG9065	Full Name: receptor (calcitonin) activity modifying protein 2	
	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean RAMP2 concentration was determined to be 2.1 ng/mL in mouse brain tissue extract based on a 1.0 mg/mL extract load, 408.1 ng/mL in mouse lung tissue extract based on a 1.0 mg/mL extract load, 1.1 ng/mL in mouse retina tissue extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03316-1, MOUSE RAMP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88020-1-PBS. Detection antibody: 88020-3-PBS. Standard: Eg9065. Range: 0.156-20 ng/mL.