

For Research Use Only

CFAP45 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88026-3-PBS



Basic Information

Catalog Number: 88026-3-PBS	GenBank Accession Number: BC027874	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 25790	CloneNo.: 260397E3
Isotype: IgG	UNIPROT ID: Q9UL16	
Immunogen Catalog Number: AG16633	Full Name: coiled-coil domain containing 19	
	Calculated MW: 551 aa, 66 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CFAP45 concentration was determined to be 85.2 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, mouse spleen tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03337-1, CFAP45 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88026-3-PBS. Detection antibody: 88026-2-PBS. Standard: Ag16633. Range: 0.078-10 ng/mL.