

For Research Use Only

PSG9 Recombinant monoclonal antibody, PBS Only

Catalog Number:88037-1-PBS



Basic Information

Catalog Number:

88037-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG8782

GenBank Accession Number:

NM_002784.5

GeneID (NCBI):

5678

UNIPROT ID:

Q00887-1

Full Name:

pregnancy specific beta-1-glycoprotein 9

Calculated MW:

48 kDa

Observed MW:

45-70 kDa

Purification Method:

Protein A purification

CloneNo.:

260068A8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88037-1-RR (PSG9 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88037-1-PBS in a different storage buffer formulation.