

For Research Use Only

PSG9 Recombinant monoclonal antibody

Catalog Number:88037-1-RR



Basic Information

Catalog Number:

88037-1-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG8782

GenBank Accession Number:

NM_002784.5

GeneID (NCBI):

5678

UNIPROT ID:

Q00887-1

Full Name:

pregnancy specific beta-1-glycoprotein 9

Calculated MW:

48 kDa

Observed MW:

45-70 kDa

Purification Method:

Protein A purification

CloneNo.:

260068A8

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human placenta tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88037-1-RR (PSG9 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.