

For Research Use Only

RNF41 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88043-3-PBS



Basic Information

Catalog Number:

88043-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11076

GenBank Accession Number:

BC032637

GeneID (NCBI):

10193

UNIPROT ID:

Q9H4P4

Full Name:

ring finger protein 41

Calculated MW:

317 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

253576D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean RNF41 concentration was determined to be 6.5 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, 9.4 ng/mL in mouse brain tissue extract based on a 6.3 mg/mL extract load, mouse eye tissue extract were served as a negative control.



Cytometric bead array standard curve of MP03321-1, RNF41 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88043-3-PBS. Detection antibody: 88043-2-PBS. Standard: Ag11076. Range: 0.625-80 ng/mL.