

For Research Use Only

SUB1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88052-2-PBS



Basic Information

Catalog Number:

88052-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1390

GenBank Accession Number:

BC009610

GeneID (NCBI):

10923

UNIPROT ID:

P53999

Full Name:

SUB1 homolog (S. cerevisiae)

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

260206E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean SUB1 concentration was determined to be 1,995 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 2,108 ng/mL in HL-60 cell extract based on a 1.2 mg/mL extract load, 254 ng/mL in Ramos cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP03355-1, SUB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88052-3-PBS. Detection antibody: 88052-2-PBS. Standard: Ag1390. Range: 0.781-100 ng/mL.