

For Research Use Only

RNASE2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88054-1-PBS



Basic Information

Catalog Number: 88054-1-PBS	GenBank Accession Number: BC093678	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6036	CloneNo.: 260304C6
Isotype: IgG	UNIPROT ID: P10153	
Immunogen Catalog Number: AG28325	Full Name: ribonuclease, RNase A family, 2 (liver, eosinophil-derived neurotoxin)	
	Calculated MW: 161 aa, 18 kDa	
	Observed MW: 18 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human urine sample were subjected to SDS PAGE followed by western blot with 88054-1-RR (RNASE2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88054-1-PBS in a different storage buffer formulation.