

For Research Use Only

KISS1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88055-3-PBS



Basic Information

Catalog Number:

88055-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27198

GenBank Accession Number:

BC022819

GeneID (NCBI):

3814

UNIPROT ID:

Q15726

Full Name:

KiSS-1 metastasis-suppressor

Calculated MW:

15 kDa

Purification Method:

Protein A purification

CloneNo.:

260303E12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean KISS1 concentration was determined to be 6.27 ng/mL in human placenta extract based on a 3.80 mg/mL extract load.



Sandwich ELISA standard curve of MP03401-1, Human KISS1 Recombinant Matched Antibody Pair - PBS only. 88055-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag27198. 88055-2-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL.