

For Research Use Only

GRB2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88060-2-PBS



Basic Information

Catalog Number: 88060-2-PBS	GenBank Accession Number: BC000631	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2885	CloneNo.: 260129D6
Isotype: IgG	UNIPROT ID: P62993	
Immunogen Catalog Number: AG0395	Full Name: growth factor receptor-bound protein 2	
	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GRB2 concentration was determined to be 16.3 ng/mL in HEK-293 cell extract based on a 1.5 mg/mL extract load, 8.2 ng/mL in THP-1 cell extract based on a 1.5 mg/mL extract load, 0.6 ng/mL in HaCaT cell extract based on a 1.4 mg/mL extract load.



Cytometric bead array standard curve of MP03345-1, GRB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88060-1-PBS. Detection antibody: 88060-2-PBS. Standard: Ag0395. Range: 0.312-40 ng/mL.