

For Research Use Only

# CSTF2T Recombinant monoclonal antibody

Catalog Number: 88063-1-RR



## Basic Information

Catalog Number:

88063-1-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5411

GenBank Accession Number:

BC028239

GeneID (NCBI):

23283

UNIPROT ID:

Q9H0L4

Full Name:

cleavage stimulation factor, 3' pre-RNA, subunit 2, 64kDa, tau variant

Calculated MW:

64 kDa

Observed MW:

~70 kDa

Purification Method:

Protein A purification

CloneNo.:

260142B6

Recommended Dilutions:

WB: 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB: HeLa cells, rat brain tissue, HepG2 cells, Jurkat cells, PC-3 cells, A549 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88063-1-RR (CSTF2T antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



rat brain tissue were subjected to SDS PAGE followed by western blot with 88063-1-RR (CSTF2T antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.