

For Research Use Only

HHLA2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88064-3-PBS



Basic Information

Catalog Number: 88064-3-PBS	GenBank Accession Number: NM_007072	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 11148	CloneNo.: 260180E5
Isotype: IgG	UNIPROT ID: Q9UM44-1	
Immunogen Catalog Number: EG5055	Full Name: HERV-H LTR-associating 2	
	Calculated MW: 414 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HHLA2 concentration was determined to be 86.2 ng/mL in A-204 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03313-1, HHLA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88064-3-PBS. Detection antibody: 88064-2-PBS. Standard: Eg5055. Range: 0.156-20 ng/mL.