

p120 Catenin Recombinant monoclonal antibody

Catalog Number: 88083-1-RR

Basic Information

Catalog Number: 88083-1-RR	GenBank Accession Number: BC010501	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 1500	CloneNo.: 260247D5
Isotype: IgG	UNIPROT ID: O60716	Recommended Dilutions: WB: 1:5000-1:50000 IHC: 1:500-1:2000
Immunogen Catalog Number: AG2824	Full Name: catenin (cadherin-associated protein), delta 1	
	Calculated MW: 948 aa, 105 kDa	
	Observed MW: 75-120 kDa	

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: HT-29 cells, HeLa cells, HEK-293 cells, MDA-MB-231 cells, NIH/3T3 cells, C6 cells

IHC: human colon tissue,

Background Information

Catenins were discovered as proteins that are linked to the cytoplasmic domain of transmembrane cadherins (PMID: 9653641). p120 catenin, also called p120 ctn or catenin delta-1 (CTNND1), regulates cell-cell adhesion through its interaction with the cytoplasmic tail of classical and type II cadherins. p120 catenin is a tyrosine kinase substrate implicated in cell transformation by SRC, as well as in ligand-induced receptor signaling through the EGF receptor, the PDGF receptor, and the CSF1 receptor. Different expression patterns of p120 catenin in lobular and ductal carcinomas of the breast have been reported: membrane stain for ductal carcinoma and cytoplasmic stain for lobular carcinoma (PMID: 24966968). Different isoforms of p120 catenin are variably expressed in different tissues as a result of alternative splicing and the use of multiple translation initiation codons (PMID: 19150613).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88083-1-RR (CTNND1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. Low expression: HL-60 cells.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 88083-1-RR (p120 Catenin antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 88083-1-RR (p120 Catenin antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).