

For Research Use Only

TRP2/DCT Recombinant monoclonal antibody, PBS Only

Catalog Number: 88100-1-PBS



Basic Information

Catalog Number: 88100-1-PBS	GenBank Accession Number: BC028311	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 1638	CloneNo.: 260396B12
Isotype: IgG	UNIPROT ID: P40126	
Immunogen Catalog Number: AG26165	Full Name: dopachrome tautomerase (dopachrome delta-isomerase, tyrosine-related protein 2)	
	Calculated MW: 519 aa, 59 kDa	
	Observed MW: 59-63 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, rat, monkey

Background Information

DCT (L-dopachrome Delta-isomerase) is also named as TYRP2, TRP2 and belongs to the tyrosinase family. Human TYRP2 protein has 83% identity and 90% similarity to the mouse sequence and has all the structural characteristics of the tyrosinase protein family, including a signal peptide, 15 conserved cysteine residues, 2 copper-binding domains, and a C-terminal membrane-spanning region (PMID:8206391). TYRP2 is involved in neural progenitor cell proliferation and has been implicated in drug and radiation resistance in malignant melanoma. Mutations in DCT cause oculocutaneous albinism type VIII (OCA8), and TYRP2 serves as both an early melanoblast marker and a melanoma-associated antigen with clinical relevance as a potential biomarker in immune checkpoint therapy.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88100-1-RR (DCT antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88100-1-PBS in a different storage buffer formulation.