

For Research Use Only

PLEKHA1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88114-1-PBS



Basic Information

Catalog Number: 88114-1-PBS	GenBank Accession Number: BC001136	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 59338	CloneNo.: 260319B3
Isotype: IgG	UNIPROT ID: Q9HB21	
Immunogen Catalog Number: AG0318	Full Name: pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 1	
	Calculated MW: 46 kDa	
	Observed MW: 46 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Pleckstrin homology (PH) domain is commonly found in eukaryotic signaling proteins and possesses multiple functions including the abilities to bind inositol phosphates and various proteins. The tandem PH domain containing protein-1 (TAPP1) or PH domain containing-family A (phosphoinositide binding specific) member 1 (PLEKHA1), interacts strongly and specifically with phosphatidylinositol 3,4-trisphosphate [PtdIns(3,4)P(2)], which is one of the immediate breakdown products of PtdIns(3,4,5) P(3) and functions as a signalling molecule in insulin- and growth-factor-stimulated pathways. TAPP1 is also associated with the protein-tyrosine-phosphatase-like protein-1 (PTPL1 also known as FAP-1) and maintains PTPL1 in cytoplasm. By binding to PtdIns(3,4) P(2) and PTPL1, TAPP1 may regulate the membrane localization of PTPL1.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88114-1-RR (PLEKHA1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88114-1-PBS in a different storage buffer formulation.