

For Research Use Only

MFI2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88127-1-PBS



Basic Information

Catalog Number:

88127-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG9197

GenBank Accession Number:

NM_005929.6

GeneID (NCBI):

4241

UNIPROT ID:

P08582-1

Full Name:

antigen p97 (melanoma associated)
identified by monoclonal antibodies
133.2 and 96.5

Calculated MW:

80 kDa

Observed MW:

82 kDa

Purification Method:

Protein A purification

CloneNo.:

260392D5

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

MFI2 (Melanotransferrin/Mtf/TRMF), also known as p97 or CD228, is an iron-binding glycoprotein belonging to the transferrin family. Alternative RNA splicing yields two forms of different length transcripts, one (the longer) bound to the membrane (mMFI2) and the other (the shorter) secreted to the peripheral environment (sMFI2) such as blood and cell supernatant. MFI2 is upregulated in tumor tissues and has been identified as a serological marker of colorectal cancer (PMID: 25216327).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



SK-MEL-30 cells were subjected to SDS PAGE followed by western blot with 88127-1-RR (MELTF/CD228 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. Lower expression of Hela cells and HEK-293 cells. This data was developed using the same antibody clone with 88127-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 88127-1-RR (MF12 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88127-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed A375 cells using MF12 antibody (88127-1-RR, Clone: 260392D5) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 88127-1-PBS in a different storage buffer formulation.