

For Research Use Only

ROBO4 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88139-2-PBS



Basic Information

Catalog Number: 88139-2-PBS	GenBank Accession Number: NM_019055.6	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 54538	CloneNo.: 260378E9
Isotype: IgG	UNIPROT ID: Q8WZ75-1	
Immunogen Catalog Number: EG7821	Full Name: roundabout homolog 4, magic roundabout (Drosophila)	
	Calculated MW: 107 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ROBO4 concentration was determined to be 94.5 ng/mL in Human plasma 1, 74.4 ng/mL in Human plasma 2.



Cytometric bead array standard curve of MP03363-1, ROBO4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88139-3-PBS. Detection antibody: 88139-2-PBS. Standard: Eg7821. Range: 0.312-40 ng/mL.