

For Research Use Only

PITX1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 88143-1-PBS



Basic Information

Catalog Number:

88143-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1303

GenBank Accession Number:

BC009412

GeneID (NCBI):

5307

UNIPROT ID:

P78337

Full Name:

paired-like homeodomain 1

Calculated MW:

34 kDa

Observed MW:

38-40 kDa

Purification Method:

Protein A purification

CloneNo.:

260497B7

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

PITX1 (Paired-like homeodomain 1) is a bicoid-related homeobox transcription factor that plays critical roles in embryonic development, particularly in patterning the hindlimb and regulating pituitary organogenesis. It is also known as "Backfoot" due to its involvement in establishing left-right asymmetry during limb formation. Beyond development, PITX1 has been implicated in various cancers, where it can function as either a tumor suppressor or an oncogene depending on the cellular context, and its dysregulation is associated with tumor progression and metastasis.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 88143-1-RR (PITX1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88143-1-PBS in a different storage buffer formulation.