

For Research Use Only

Mouse HA2U Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88145-3-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 88145-3-PBS | GenBank Accession Number: / | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 14960 | CloneNo.: 260473B9 |
| Isotype: IgG | UNIPROT ID: P14438 | |
| Immunogen Catalog Number: EG9227 | Full Name: histocompatibility 2, class II antigen A, alpha | |
| | Calculated MW: 25 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HA2U concentration was determined to be 183.3 ng/mL in mouse spleen tissue extract based on a 1.0 mg/mL extract load, 109.7 ng/mL in mouse thymus tissue extract based on a 1.0 mg/mL extract load, 1.8 ng/mL in mouse heart tissue extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03392-1, MOUSE HA2U Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88145-4-PBS. Detection antibody: 88145-3-PBS. Standard: Eg9227. Range: 0.156-20 ng/mL.