

For Research Use Only

FCRL5 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88157-1-PBS



Basic Information

Catalog Number: 88157-1-PBS	GenBank Accession Number: BC101067	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 83416	CloneNo.: 253566C6
Isotype: IgG	UNIPROT ID: Q96RD9	
Immunogen Catalog Number: EG1902	Full Name: Fc receptor-like 5	
	Calculated MW: 977 aa, 106 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FCRL5 concentration was determined to be 34.7 ng/mL in Human plasma 1, 46.7 ng/mL in Human plasma 2.



Cytometric bead array standard curve of MP03348-1, FCRL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88157-2-PBS. Detection antibody: 88157-1-PBS. Standard: Eg1902. Range: 0.625-80 ng/mL.