

For Research Use Only

# FCRL5 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88157-2-PBS



## Basic Information

<b>Catalog Number:</b> 88157-2-PBS	<b>GenBank Accession Number:</b> BC101067	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 83416	<b>CloneNo.:</b> 253566C2
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> Q96RD9	
<b>Immunogen Catalog Number:</b> EG1902	<b>Full Name:</b> Fc receptor-like 5	
	<b>Calculated MW:</b> 977 aa, 106 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean FCRL5 concentration was determined to be 34.7 ng/mL in Human plasma 1, 46.7 ng/mL in Human plasma 2.



Cytometric bead array standard curve of MP03348-1, FCRL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88157-2-PBS. Detection antibody: 88157-1-PBS. Standard: Eg1902. Range: 0.625-80 ng/mL.