

For Research Use Only

G6B Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88167-1-PBS



Basic Information

Catalog Number: 88167-1-PBS	GenBank Accession Number: NM_138272.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 80739	CloneNo.: 260393E3
Isotype: IgG	UNIPROT ID: O95866-1	
Immunogen Catalog Number: EG5762	Full Name: chromosome 6 open reading frame 25	
	Calculated MW: 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MPIG6B concentration was determined to be 8.8 ng/mL in HEL cell extract based on a 1.5 mg/mL extract load, THP-1 cell extract were served as a negative control.



Cytometric bead array standard curve of MP03352-1, MPIG6B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88167-2-PBS. Detection antibody: 88167-1-PBS. Standard: Eg5762. Range: 0.156-20 ng/mL