

For Research Use Only

G6B Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88167-3-PBS



Basic Information

Catalog Number:

88167-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG5762

GenBank Accession Number:

NM_138272.3

GeneID (NCBI):

80739

UNIPROT ID:

O95866-1

Full Name:

chromosome 6 open reading frame 25

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

260393E9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MPIG6B concentration was determined to be 2.2 ng/mL in K-562 cell extract based on a 1.2 mg/mL extract load, 5.8 ng/mL in HEL cell extract based on a 1.5 mg/mL extract load, 1.1 ng/mL in 786-O cell extract based on a 1.2 mg/mL extract load, 2.4 ng/mL in THP-1 cell extract based on a 1.5 mg/mL extract load and 1.1 ng/mL in RT-4 cell extract based on a 1.5 mg/mL extract load.



Sandwich ELISA standard curve of MP03352-2, Human MPIG6B Recombinant Matched Antibody Pair - PBS only. 88167-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg5762. 88167-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.