

For Research Use Only

ACD Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88168-3-PBS



Basic Information

Catalog Number: 88168-3-PBS	GenBank Accession Number: BC016904	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 65057	CloneNo.: 260314F11
Isotype: IgG	UNIPROT ID: Q96APO	
Immunogen Catalog Number: AG23021	Full Name: adrenocortical dysplasia homolog (mouse)	
	Calculated MW: 58 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ACD concentration was determined to be 7.4 ng/mL in Jurkat cell extract based on a 1.0 mg/mL extract load, 8.7 ng/mL in HeLa cell extract based on a 1.0 mg/mL extract load, 4.3 ng/mL in SK-BR-3 cell extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03378-1, ACD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88168-3-PBS. Detection antibody: 88168-2-PBS. Standard: Ag23021. Range: 0.312-40 ng/mL