

For Research Use Only

CCL5/RANTES Recombinant monoclonal antibody, PBS Only

Catalog Number: 88173-1-PBS



Basic Information

Catalog Number:

88173-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2617

GenBank Accession Number:

BC008600

GeneID (NCBI):

6352

UNIPROT ID:

P13501

Full Name:

chemokine (C-C motif) ligand 5

Calculated MW:

91 aa, 10 kDa

Observed MW:

8 kDa

Purification Method:

Protein A purification

CloneNo.:

260534G8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

RANTES (CCL5) belongs to the CC chemokine family and induces leukocyte migration by binding to specific receptors in the G protein-coupled receptor family. RANTES (CCL5) has been shown in vitro to mediate eosinophil, lymphocyte, neutrophil, and monocyte chemotaxis, and it initiates several other proinflammatory events, such as integrin activation, lipid mediator biosynthesis, and degranulation.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human peripheral blood platelets were subjected to SDS PAGE followed by western blot with 88173-1-RR (CCL5/RANTES antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88173-1-PBS in a different storage buffer formulation.



PMA, LPS and Brefeldin A-treated THP-1 cells were subjected to SDS PAGE followed by western blot with 88173-1-RR (CCL5/RANTES antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 88173-1-PBS in a different storage buffer formulation.