

For Research Use Only

# S100P Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88256-1-PBS



## Basic Information

Catalog Number:

88256-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2379

GenBank Accession Number:

BC006819

GeneID (NCBI):

6286

UNIPROT ID:

P25815

Full Name:

S100 calcium binding protein P

Calculated MW:

95 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

260786G1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean S100P concentration was determined to be 464.9 ng/mL in BxPC-3 cell extract based on a 1.2 mg/mL extract load, SH-SY5Y cell extract were served as a negative control.



Cytometric bead array standard curve of MP03518-1, S100P Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88256-1-PBS. Detection antibody: 88256-3-PBS. Standard: Ag2379. Range: 0.312-40 ng/mL