

For Research Use Only

S100P Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88256-3-PBS



Basic Information

Catalog Number: 88256-3-PBS	GenBank Accession Number: BC006819	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6286	CloneNo.: 260786D7
Isotype: IgG	UNIPROT ID: P25815	
Immunogen Catalog Number: AG2379	Full Name: S100 calcium binding protein P	
	Calculated MW: 95 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean S100P concentration was determined to be 464.9 ng/mL in BxPC-3 cell extract based on a 1.2 mg/mL extract load, SH-SY5Y cell extract were served as a negative control.



Cytometric bead array standard curve of MP03518-1, S100P Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88256-1-PBS. Detection antibody: 88256-3-PBS. Standard: Ag2379. Range: 0.312-40 ng/mL