

For Research Use Only

CYP4F11 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 88260-3-PBS



Basic Information

Catalog Number:

88260-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2906

GenBank Accession Number:

BC016853

GeneID (NCBI):

57834

UNIPROT ID:

Q9HBI6

Full Name:

cytochrome P450, family 4, subfamily F, polypeptide 11

Calculated MW:

524 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

260726D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CYP4F11 concentration was determined to be 3.0 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load, Ramos cell extract were served as a negative control.



Cytometric bead array standard curve of MP03443-1, CYP4F11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88260-3-PBS. Detection antibody: 88260-2-PBS. Standard: Ag2906. Range: 0.312-40 ng/mL.