

For Research Use Only

NOC2L Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 88267-2-PBS



Basic Information

Catalog Number:

88267-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26952

GenBank Accession Number:

BC003555

GeneID (NCBI):

26155

UNIPROT ID:

Q9Y3T9

Full Name:

nucleolar complex associated 2 homolog (S. cerevisiae)

Calculated MW:

749 aa, 85 kDa

Purification Method:

Protein A purification

CloneNo.:

260756E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NOC2L concentration was determined to be 26.5 ng/mL in Raji cell extract based on a 1.0 mg/mL extract load, 26.8 ng/mL in HT-1080 cell extract based on a 1.0 mg/mL extract load, 20.3 ng/mL in U2OS cell extract based on a 1.0 mg/mL extract load.



Cytometric bead array standard curve of MP03442-1, NOC2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 88267-1-PBS. Detection antibody: 88267-2-PBS. Standard: Ag26952. Range: 0.156-20 ng/mL.